

SAFCONET TABS®

Micro-organisms concentrate for biological parts-washer

DESCRIPTION

SAFCONET TABS® are composed of a selection of micro-organisms that are specifically suitable for the bioremediation of grease and oil.

INSTRUCTIONS FOR USE

When installing your degreasing fountain, put a **SAFCONET TABS®** into the fountain's reservoir to start the process.

Add a tab every 4 to 6 weeks, depending on conditions of use of your fountain.

ADVANTAGES

- Non-irritant for skin & eyes.
- Non-caustic.
- No labelling.
- Neutralizes odors.
- Reduces sludge.

TECHNICAL INFORMATION

TECHNICAL CHARACTERISTICS	UNITS	SAFCONET TABS®
Aspect & Color	-	Solid Blue
Temperature of use	°C	20-45
Bacteria concentration (SM9215C)	UFC/tab	20 billion
Weight	g	28.4
Solubility in water	%	100

These characteristics are given for information and can be modified over the time.

HEALTH & SAFETY

- Selection of microorganisms not belonging to the list of biological pathogens in Groups 2, 3 and 4, set by the French decree of 18 July 1994 and amended by decree of 17 April 1997 and June 30, 1998
- Non-classified as irritants for eyes and skin in accordance with criteria set forth by Directive 67/548/EEC
- No labelling
- Store in a cool, dry place away from light and heat

SERVICES AND EQUIPMENT

In addition to its product ranges, MotulTech can provide tools and services for the maintenance and in-use monitoring of your lubricants.